

# COMPONENTS CATALOGUE


FOR I.V. GIVING SETS,
BLOOD GIVING SETS,
AND OTHER MEDICAL
DEVICES


### INDEX


I.V. SET SMALL DRIP CHAMBERS	PG.
I.V. SET BIG DRIP CHAMBERS	3
In-Line Drip Chambers	5
BLOOD GIVING SET DRIP CHAMBERS	5
SPIKES AND STEPPED CONNECTOR	6
CAPS FOR MALE CONNETCTOR AND MALE CONNECTOR	6
NEEDLESS CONNECTORS	<b>7</b>
CONNECTORS WITH NON-RETURN VALVE	8
TUBING CONNECTORS	8
CONNECTORS FOR PUMP	9
INJECTION PORTS	10
CLAMPS (ROLLER AND SLIDE)	10
PATIENT SIDE TERMINALS	1 1
Butterflies	11
FEMALE CONNECTOR AND FEMALE CONNECTOR CAPS	12
Burettes	13
3-WAY STOPCOCKS	13
Manifolds	14
Tubings	15
FLOW REGULATORS	16
COMPONENTS FOR ENTERAL NUTRITION	16
WITHDRAWAL CANNULAE	17
SCISSOR CLAMPS AND FORCEPS	17
BREAK-OFF TIP CONNECTORS / VIAL STOPPERS	17
IV In-LINE FILTERS	18










































#### IN-LINE IV ADULT SQUARE FILTERS WITH 10 CM OF FILTRATION AREA 5 MICRON CLEAR FILTER FIT TO (PVC) (PVC) TUBE 3.0 x 4.1 MM 3.0 x 4.1 MM FOR PHOTOSENSITIVE DRUGS FILTRATION AREA: 10 CM2 FILTRATION AREA: 10 CM2 HOUSING: CLEAR MODIFIED ACRYLIC 5 micron HOUSING: AMBER/OPAQUE micron MODIFIED ACRYLIC SOLUTION FILTER: HYDROPHILIC SOLUTION FILTER: HYDROPHILIC PES VENT FILTER: HYDROPHOBIC VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (0.03 MICRON) PTFE MEMBRANE (0.03 MICRON) MINIMUN FLOW RATE: 95ML/MIN AT 80 CM WATER HEAD MINIMUN FLOW RATE: 95ML/MIN PRESSURE AT 80 CM WATER HEAD PRESSURE CODE: 103000 CARTON QUANTITY: 3.600 CODE: 1C30000P CARTON QUANTITY: 3.600 1.2 MICRON FILTER FIT TO TUBE 1.2 MICRON CLEAR FILTER FIT (PVC) PVC TO TUBE 3.0 x 4.1 MM 3.0 x 4.1 MM FOR PHOTOSENSITIVE DRUGS FILTRATION AREA: 10 CM FILTRATION AREA: 10 CM<sup>2</sup> HOUSING: CLEAR MODIFIED micron micror ACRYLIC HOUSING: AMBER/OPAQUE MODIFIED ACRYLIC SOLUTION FILTER: HYDROPHILIC 1.2 7 PES SOLUTION FILTER: HYDROPHILIC PES VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (0.03 MICRON) VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (0.03 MICRON) MINIMUN FLOW RATE: 52 ML/MIN MINIMUN FLOW RATE: 52 ML/MIN AT 80 CM WATER HEAD PRESSURE AT 80 CM WATER HEAD PRESSURE CODE: ICSOOS CARTON QUANTITY: 3.600 CARTON QUANTITY: 3.600 CODE: ICSCOZOR D.2 MICRON CLEAR FILTER FIT D.2 MICRON FILTER FIT TO TUBE TO TUBE 3.0 X 4.1 MM 3.0 x 4.1 MM FOR PHOTOSENSITIVE DRUGS FILTRATION AREA: 10 CM FILTRATION AREA: 10 CM<sup>2</sup> HOUSING: CLEAR MODIFIED micron micron HOUSING: AMBER/OPAQUE MODIFIED ACRYLIC SOLUTION FILTER: HYDROPHILIC PES SOLUTION FILTER: HYDROPHILIC PES VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (0.03 MICRON) VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (0.03 MICRON) MINIMUN FLOW RATE: 20 ML/MIN AT 80 CM WATER HEAD MINIMUN FLOW RATE: 20 ML/MIN PRESSURE AT 80 CM WATER HEAD PRESSURE CODE: 103001 CARTON QUANTITY: 3.600 CODE: 1030010P CARTON QUANTITY: 3.600 OTHER TECHNICAL SPECIFICATION: -DIMENSION: WIDTH 30.0 MM, DEPTH 6.5 MM, LENGTH 60 MM (INCLUDING TUBE CONNECTORS), 40 MM (FILTER BODY) -WEIGHT: 5.2 GR -EFFECTIVE FILTRATION: 99.9999% -REF. STANDARD: ISO 8536-4 -MINIMUN WATER BUBBLE POINT: 3.7-4.8 BAR MAXIMUN DEFRATING TEMPERATURE: 5-40°C -PRIMING VOLUME >2.4 ML



#### IN-LINE IV ADULT SQUARE FILTERS WITH 5 CM OF FILTRATION AREA 5 MICRON CLEAR FILTER FIT TO 5 MICRON FILTER FIT TO TUBE (PVC) (PVC) TUBE 3.0 x 4.1 MM 3.0 x 4.1 MM FOR PHOTOSENSITIVE DRUGS FILTRATION AREA: 5 CM2 FILTRATION AREA: 5 CM2 HOUSING: CLEAR PP 5 micron HOUSING: AMBER/OPAQUE PP micron SOLUTION FILTER: HYDROPHILIC SOLUTION FILTER: HYDROPHILIC PES VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (0.5 MICRON) VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (0.5MICRON) MINIMUN FLOW RATE: 50ML/MIN AT 80 CM WATER HEAD MINIMUN FLOW RATE: 50ML/MIN PRESSURE AT 80 CM WATER HEAD PRESSURE CODE: 103202 CARTON QUANTITY: 4.000 CODE: IC32020P CARTON QUANTITY: 4.000 1.2 MICRON FILTER FIT TO TUBE 1.2 MICRON CLEAR FILTER FIT (PVC) PVC TO TUBE 3.0 x 4.1 MM 3.0 x 4.1 MM FOR PHOTOSENSITIVE DRUGS FILTRATION AREA: 5 CM FILTRATION AREA: 5 CM2 HOUSING: CLEAR PP micron micror HOUSING: AMBER/OPAQUE PP SOLUTION FILTER: HYDROPHILIC PES SOLUTION FILTER: HYDROPHILIC 7 PES VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (D.5 MICRON) VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (O.5 MICRON) MINIMUN FLOW RATE: 50 ML/MIN AT 80 CM WATER HEAD MINIMUN FLOW RATE: 50 ML/MIN PRESSURE AT 80 CM WATER HEAD PRESSURE CODE: 103201 CARTON GUANTITY: 4.000 CODE: 1032010P CARTON GUANTITY: 4.000 D.2 MICRON CLEAR FILTER FIT D.2 MICRON FILTER FIT TO TUBE TO TUBE 3.0 X 4.1 MM 3.0 x 4.1 MM FOR PHOTOSENSITIVE DRUGS FILTRATION AREA: 5 CM FILTRATION AREA: 5 CM2 HOUSING: CLEAR PP micron micron HOUSING: AMBER/OPAQUE PP SOLUTION FILTER: HYDROPHILIC PES SOLUTION FILTER: HYDROPHILIC 0.2 PES VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (O.5 MICRON) VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (O.5 MICRON) MINIMUN FLOW RATE: 10 ML/MIN AT 80 CM WATER HEAD MINIMUN FLOW RATE: 10 ML/MIN PRESSURE AT 80 CM WATER HEAD PRESSURE CODE: 103200 CARTON QUANTITY: 4.000 CODE: 1032000P CARTON QUANTITY: 4.000 OTHER TECHNICAL SPECIFICATION: -DIMENSION: WIDTH 29.0 MM, DEPTH 8.3 MM, LENGTH 40 MM (FILTER BODY) -WEIGHT: 5.6 GR -EFFECTIVE FILTRATION: 98% -REF. STANDARD: ISO 8536-4 MAXIMUN OPERATING TEMPERATURE: 5-40°C

#### IN-LINE IV ADULT/PAEDIATRIC ROUND FILTERS

5 MICRON CLEAR FILTER FIT TO TUBE 3.0 x 4.1 MM

FILTRATION AREA: 4.2 CM2

HOUSING: CLEAR PP

SOLUTION FILTER: HYDROPHILIC PES

VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (0.5 MICRON)

MINIMUN FLOW RATE: 50ML/MIN AT 80 CM WATER HEAD PRESSURE



3.0 x 4.1 MM FOR PHOTOSENSITIVE DRUGS


FILTRATION AREA: 4.2 CM2

HOUSING: AMBER/OPAQUE PP

SOLUTION FILTER: HYDROPHILIC PES

VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (O.5MICRON)

MINIMUN FLOW RATE: 50ML/MIN AT 80 CM WATER HEAD PRESSURE



CODE: 103205

CARTON QUANTITY: 5.000

CODE: 1032050P

CARTON QUANTITY: 5.000

1.2 MICRON CLEAR FILTER FIT TO TUBE 3.0 x 4.1 MM


FILTRATION AREA: 4.2 CM

HOUSING: CLEAR PP

SOLUTION FILTER: HYDROPHILIC PES

VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (D.5 MICRON)

MINIMUN FLOW RATE: 50 ML/MIN AT 80 CM WATER HEAD PRESSURE



1.2 MICRON FILTER FIT TO TUBE 3.0 x 4.1 MM FOR PHOTOSENSITIVE DRUGS

FILTRATION AREA: 4.2 CM2

HOUSING: AMBER/OPAQUE PP

SOLUTION FILTER: HYDROPHILIC PES

VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (D.5 MICRON)

MINIMUN FLOW RATE: 50 ML/MIN AT 80 CM WATER HEAD PRESSURE



CODE: 103204

CARTON QUANTITY: 5.000

CODE: 1032040P

CARTON QUANTITY: 5.000

O.2 MICRON CLEAR FILTER FIT TO TUBE 3.0 x 4.1 MM


FILTRATION AREA: 4.2 CM2

HOUSING: CLEAR PP

SOLUTION FILTER: HYDROPHILIC PES

VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (O.5 MICRON)

MINIMUN FLOW RATE: 10 ML/MIN AT 80 CM WATER HEAD PRESSURE



O.2 MICRON FILTER FIT TO TUBE 3.0 x 4.1 MM FOR PHOTOSENSITIVE DRUGS

FILTRATION AREA: 4.2 CM<sup>2</sup>

HOUSING: AMBER/OPAQUE PP

SOLUTION FILTER: HYDROPHILIC PES

VENT FILTER: HYDROPHOBIC PTFE MEMBRANE (O.5 MICRON)

MINIMUN FLOW RATE: 10 ML/MIN AT 80 CM WATER HEAD PRESSURE



CODE: 103203

CARTON QUANTITY: 5.000

CODE: 1032030P

CARTON QUANTITY: 5.000

OTHER TECHNICAL SPECIFICATION:

-DIMENSION: DIAMETER 31 MM (FILTER BODY)

-WEIGHT: 3.8 GR

-EFFECTIVE FILTRATION: 98%

-REF. STANDARD: ISO 8536-4

MAXIMUN OPERATING TEMPERATURE: 5-40°C



















## PR o ACTIVE



LIMITS OF RESPONSIBILITY - It is your responsibility to validate the proposed material. We guarantee only the dimensional characteristics of the part and under no circumstances does our liability extend to the end product in wich it is found. In the case of dimensional non-compliance, we shall take back the parts in question and exchange them. We shall not be held responsible for the later use of the supplied pieces. We cannot guarantee the behavior of the piece with ragard to: gluing, sealing, welding (ultrasound or high frequency) nor for the different means of sterilization used: gas gamma ray, beta ray, accellerated electrons, dry or humyd steams, etc. The same applies to the different chemical, biochemical, substances or biological liquids in contact with the part. It is your responsability to check compliance of the material and the dimension of the part with the processes it shall undergo and the required purpose.

#### PRO-ACTIVE S.r.I.

ADDRESS: VIALE DEL COMMERCIO, 12 - 45100 ROVIGO - ITALY

**TELEPHONE:** +39.0425.491087 - +39.0425.1885138

FAX: +39.0425.490101

#### **E-MAIL**

**GENERAL INFORMATIONS:** info@proactivesrl.it

COMMERCIAL AND ADMINISTRATION INFORMATIONS: barbara.bozzolan@proactivesrl.it

TECHNICAL INFORMATIONS AND REGULATORY AFFAIRS: fabio.mazzocchi@proactivesrl.it

#### **COMPANY DATA**

CAP. STOCK 90.000,00 €

**VAT:** IT01248310292

BANK CODE: Cassa di Risparmio del Veneto - IBAN: IT 33 G 06225 12203 100000001146 - BIC: IBSP IT 2P